

VINYL TECHNICAL CHARACTERISTICS

- Antimicrobial protection
- Antibacterial protection
- Antimycotic protection
- Abrasion protection
 - Martindale test (dry) > 51.200 rubs
 - Martindale test (wet) > 12.800 rubs
 - Wyzenbeek test > rubs
- Resistance to spitting and perspiration
- Urine resistance
- Blood resistance
- Stain resistance surface
- Flame retardant M-2

CARE AND CLEANING GUIDE

- The vinyl should be cleaned periodically in order to maintain its appearance and prevent build-up of dirt and contaminants.
- Any stain, spills or soiling should be cleaned up promptly to prevent the possibility of permanent staining.
- Use soft soapy solutions or special cleaning products for vinyl fabrics to remove stains on the surface of the material. Remove only with a damp white cloth.
- Lacquers, strong cleaners, detergents, xylene-based solutions, acetone or ketone (MEK) cause immediate damage and contribute to the deterioration of the material. The use of such cleaners is at the owner's risk.
- Do not use abrasives, solvents or silicone compounds. Also avoid cleaning with stiff bristled brushed or metal sponges, as they can affect the material properties.

DAILY CLEANING AND MAINTENANCE

1. Clean the entire surface using of 1:9 liquid soap (neutral PH) and water. Do not use soap or cleaning solutions that contains alcohol, ketones, xylene, acetates or solvents (mineral/white spirits)
2. Completely remove excess cleaner with a white, clean, damp cloth.
3. Dry the surface.

CLEANING METHODS

Method 1: Coffee, juice, red wine, sauces, chocolate, grease, crayon and suntan lotion.

- Remove excess spill with a damp, white cloth or soft brush.
- Clean with a 1:9 mix of liquid neutral Ph soap and water.
- Rinse with clean water and dry.

Method 2: Make-up

- Use a straight application of neutral Ph soap.
- Wipe with a clean white cloth or soft brush, rinse with clean water and dry

Method 3: Blood, urine, feces

- Use a 1:1 mix of Isopropyl alcohol and water.
- Rinse with clean water and dry.

Method 4: All other tough stains

- Use a 0.5% solution of household bleach and water (1:9 mix of 5% concentrated household bleach and water).
- Rinse thoroughly with clean water and dry

Note: Do not use steel bristle brushes, abrasive scouring pads or metal, these cause permanent damage to the fabric

Method 2 and 4 are not recommended for every day cleaning!

This information is not a guarantee and does not relieve the user from the responsibility of the proper and safe use of the product and all cleaning agents. The use of certain agents can be harmful to the surface appearance and lifespan of vinyl. INMOCLINC will not assume responsibility resulting from the use of such cleaning agents to the vinyl.

Certain clothing and accessory dyes (such as those used on denim jeans) may migrate to lighter colours. This phenomenon is increased by humidity and temperature and is irreversible.

INMOCLINC is not the manufacturer of this vinyl. No responsibility is taken for the correctness of this information. All information is according to the manufacturers indications.

Colors

